

GIROFLEX ADAPT



REHA



giroflex
designed to work

Developed for coxarthrosis patients

Elevated seat position and seat inclination towards the front

Excellent sitting comfort

Available with climate-functional Giroflex “Grano” cover fabric

Available in artificial leather and various colours for clinical environments

Development in cooperation with doctors and patients

GIROFLEX ADAPT REHA-COXARTHROSIS ARMCHAIR

The helpful ones. Many people suffer from arthritis in the hip area. The impairment is extremely painful. Above all, sitting after a hip joint surgery is often very painful for the person who had the operation. The giroflex 33 has been specially developed for the needs of coxarthrosis patients. The armchair makes it possible for the affected person to sit without complaints due to an elevated sitting position and a seat angle that is specially inclined towards the front. The definitive seat height also contributes to the greatest possible sitting comfort under special condi-

tions. This helps to relieve the pelvis by leaning forward. Upon request, the giroflex Reha-Coxarthrosis Armchair can be equipped with castors. This gives it mobility and allows it to be moved more easily. The giroflex 33 is equally suited for use in therapeutic facilities and private environments. In its basic version, the Reha-Coxarthrosis Armchair already has a chromium-plated base frame. As a result, it is easier to clean, looks attractive and has high life expectancy.



Optionally available with crutch holder



33-7008-RC

- Artificial leather cover
- Sled-base frame, chromed



33-7009-RC

- Artificial leather cover
- Sled-base frame, chromed
- Castors 50 mm



33-7009-RC

- Artificial leather cover
- Crutch holder
- Sled-base frame, chromed
- Castors 50 mm

GIROFLEX ADAPT REHA-COXARTHROSIS ARMCHAIR



33-7008-RC



33-7009-RC

Dimensions in cm		
Height after use	53	53
Seat tilt forwards	-3°	-3°
Total height	93	93
Chairback height	38	38
Armrest height	78	78
Seat width	42	42
Chairback width	44	44
Total width	64	64
Seat depth	43	43

giroflex adapt REHA-Coxarthrosis Armchair: Cover materials made of artificial leather skai@Palma or fabric "Strong". Optionally available with crutch holders, gliders or 50 mm castors. These options are always just reserved for the specific models.

Switzerland

Stoll Giroflex AG
CH-5322 Koblenz
Tel.: +41 (0) 56 267 91 11
info@giroflex.ch

Germany

Giroflex GmbH
D-78187 Geisingen
Tel.: +49 (0) 7704 9279 0
info@giroflex.de

Belgium

Giroflex SA
B-1200 Brussel
Tel.: +32 (0) 2 761 20 20
info@giroflex.be

Netherlands

Giroflex Netherlands BV
NL-1184 VW Ouderkerk
aan de Amstel
Tel.: +31 (0) 20 47 22 555
info@giroflex.nl

France

Giroflex France
Tel.: +33 (0) 800 91 72 71
info@giroflex.be

www.giroflex.com

giroflex
designed to work